



atro
Flexible Implant



Flexible Neck Implant

Alfa gate Patro implant are specifically designed for narrow ridges and tight spaces.

The insertion of Patro implant is quick and simple one stage procedure, these versatile implant can be used for single or multiple unite restoration with immediate loading in both upper and lower jaws.

The flexible neck specially designed for immediate abutment angulation adjustment.



DETERMINED TO LEAD 



www.alfagate.com


Cycon Catalog

AlfaGate
IMPLANTS YOU CAN TRUST

atro
Flexible Implant



Patro offers the perfect solution:
for a unique one-piece implant with its flexible neck design that reaches up to 25° angulation.




The implant's versatility:
offers immediate loading in multiple-unit restorations



Conveniently it can be:
placed in the tightest and narrowest ridges.




Diameters




Ø mm	L mm	6	8	10	11.5	13	16
3.0				✓	✓	✓	✓
3.3			✓	✓	✓	✓	✓
3.75			✓	✓	✓	✓	✓
4.2		✓	✓	✓	✓	✓	✓
4.7		✓	✓	✓	✓	✓	✓

A wide range of diameters from 3.0mm up to the special 4.7mm.

Body Characteristics



Flexible Neck



Sharp and Deep Threads

Material:

Titanium Alloy Ti6Al4V ELI

Item includes:

Cover screw and implant carrier

Connection:

One-Piece implant with an integrated abutment

Surface:

Sand blasting and Acid-Etching